

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CXCR5 Polyclonal Antibody

多克隆抗体

产品货号: PA06536

第三版

描述: Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes. May have a regulatory function in Burkitt lymphoma (BL) lymphomagenesis and/or B-cell differentiation.

抗原: Synthesized peptide derived from the N-terminal region of human Blr1.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expression in mature B-cells and Burkitt lymphoma cells.