

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD240d Polyclonal Antibody

多克隆抗体

产品货号: PA06545

第三版

描述:The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene, which encodes the RhD protein, and a second gene that encodes both the RhC and RhE antigens on a single polypeptide. The two genes, and a third unrelated gene, are found in a cluster on chromosome 1. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. Multiple transcript variants encoding different isoforms have been found for this gene. RHD (Rh Blood Group D Antigen) is a Protein Coding gene. Diseases associated with RHD include Anemia, Hemolytic, Rh-Null, Regulator Type and Skeletal Muscle Cancer. GO annotations related to this gene include ammonium transmembrane transporter activity. An important paralog of this gene is RHCE.

抗原:Synthesized peptide derived from the Internal region of human CD240d.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Restricted to tissues or cell lines expressing erythroid characters.