

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CD232 Polyclonal Antibody

多克隆抗体

产品货号: PA06546

第三版

描述:This gene encodes a member of the plexin family. Plexins are transmembrane receptors for semaphorins, a large family of proteins that regulate axon guidance, cell motility and migration, and the immune response. The encoded protein and its ligand regulate melanocyte adhesion, and viral semaphorins may modulate the immune response by binding to this receptor. The encoded protein may be a tumor suppressor protein for melanoma. Alternatively spliced transcript variants have been observed for this gene. PLXNC1 (Plexin C1) is a Protein Coding gene. Diseases associated with PLXNC1 include Polycystic Kidney And Hepatic Disease. Among its related pathways are Developmental Biology and Semaphorin interactions. GO annotations related to this gene include receptor binding and receptor activity. An important paralog of this gene is PLXND1.

抗原:Synthesized peptide derived from the Internal region of human CD232.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Detected in heart, brain, lung, spleen and placenta.