

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GDF3 Polyclonal Antibody

多克隆抗体

产品货号: PA06563

第三版

描述:The protein encoded by this gene is a member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily. This group of proteins is characterized by a polybasic proteolytic processing site which is cleaved to produce a mature protein containing seven conserved cysteine residues. The members of this family are regulators of cell growth and differentiation in both embryonic and adult tissues.

抗原:Synthesized peptide derived from the C-terminal region of human GDF-3.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: