
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NRTN Polyclonal Antibody

多克隆抗体

产品货号: PA06566

第三版

描述:Neurturin is a member of the TGF-beta subfamily, TRN. This gene signals through RET and a GPI-linked coreceptor, and promotes survival of neuronal populations. A neurturin mutation has been described in a family with Hirschsprung Disease.

抗原:Synthesized peptide derived from the Internal region of human NTN.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: