

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EPGN Polyclonal Antibody

多克隆抗体

产品货号: PA06573

第三版

描述:The protein encoded by this gene is a member of the epidermal growth factor family. Members of this family are ligands for the epidermal growth factor receptor and play a role in cell survival, proliferation and migration. This protein has been reported to have high Human, Mouse mitogenic activity but low affinity for its receptor. Expression of this transcript and protein have been reported in cancer specimens of the breast, bladder, and prostate. Alternative splicing results in multiple transcript variants.

抗原:Synthesized peptide derived from the Internal region of human Epigen.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: