

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat FUNDC1 Polyclonal Antibody

多克隆抗体

产品货号: PA06575

第三版

描述: This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known. FUNDC1 (FUN14 Domain Containing 1) is a Protein Coding gene. Diseases associated with FUNDC1 include Hypoxia. Among its related pathways are Mitophagy - animal. An important paralog of this gene is FUNDC2.

抗原: Synthesized peptide derived from the N-terminal region of human FUNDC1.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed.