

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD66e Polyclonal Antibody

多克隆抗体

产品货号: PA06578

第三版

描述: Carcinoembryonic antigen-related cell adhesion molecule 5 (CEACAM5) also known as CD66e (Cluster of Differentiation 66e), is a human gene. Members of the CEACAM subfamily, including CEACAM5, belong to the CEA gene family. CD66E is a well known tumor marker and a heavily glycosylated GPI-linked cell surface molecule.

抗原: Recombinant protein of human CEACAM5

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.