

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TONSL Polyclonal Antibody

多克隆抗体

产品货号: PA06650

第三版

描述:The protein encoded by this gene is thought to be a negative regulator of NF-kappa-B mediated transcription. The encoded protein may bind NF-kappa-B complexes and trap them in the cytoplasm, preventing them from entering the nucleus and interacting with the DNA. Phosphorylation of this protein targets it for degradation by the ubiquitination pathway, which frees the NF-kappa-B complexes to enter the nucleus.

抗原:Recombinant protein of human TONSL

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart, skeletal muscle and tracheal epithelial cells.