

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat BMP3 Polyclonal Antibody

多克隆抗体

产品货号: PA06736

第三版

描述:BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation. Negatively regulates bone density. Antagonizes the ability of certain osteogenic BMPs to induce osteoprogenitor differentiation and ossification. Expressed in adult and fetal cartilage. Highly expressed in fracture tissue, particularly in osteoblasts, osteoclasts and chondroblasts.

抗原:Recombinant protein of human BMP3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in adult and fetal cartilage.