

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BTG3 Polyclonal Antibody

多克隆抗体

产品货号: PA06740

第三版

描述:The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein might play a role in neurogenesis in the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.

抗原:Recombinant protein of human BTG3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous. High expression in the ventricular zone of the developing central nervous system. High in ovary, testis, prostate, thymus and lung.