

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CNNM3 Polyclonal Antibody

多克隆抗体

产品货号: PA06813

第三版

描述:The ACDP3 gene, designated CNNM3 by the HUGO Nomenclature Committee, encodes a deduced 707-amino acid protein. Northern blot analysis revealed ubiquitous expression of CNNM3. Immunofluorescence studies showed that all 4 CNMM proteins are predominantly localized in the nucleus. Acdp3 contains 4 transmembrane domains, 2 CBS domains, a DUF21 domain, and a cNMP-binding domain. It also has large alanine- and leucine-rich regions. Northern blot analysis detected variable Acdp3 expression in all mouse tissues examined, with highest expression in brain, kidney, liver, and heart.

抗原:Recombinant protein of human CNNM3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed. Expressed at higher level in heart and spleen.