

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EDN2 Polyclonal Antibody

多克隆抗体

产品货号: PA06851

第三版

描述:Endothelin 2, also known as ET-2, is a protein that in humans is encoded by the EDN2 gene. This gene encodes a member of the endothelin protein family of secretory vasoconstrictive peptides. The preproprotein is processed to a short mature form which functions as a ligand for the endothelin receptors that initiate intracellular signaling events. This gene product is involved in a wide range of biological processes, such as hypertension and ovulation.

抗原:Recombinant protein of human EDN2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in lung, but not in placental stem villi vessels or cultured placental villi smooth muscle cells.