

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Galectin 8 Polyclonal Antibody

多克隆抗体

产品货号: PA06873

第三版

描述: This gene encodes a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions. Alternatively spliced transcript variants encoding different isoforms have been identified.

抗原: Recombinant protein of human LGALS8

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous. Selective expression by prostate carcinomas versus normal prostate and benign prostatic hypertrophy.