

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse KLF15 Polyclonal Antibody

多克隆抗体

产品货号: PA06947

第三版

描述: Kr ü ppe-like factor 15 is a protein that in humans is encoded by the KLF15 gene in the Kr ü ppe-like factor family. KLF15, a kidney-enriched Kr ü ppe-like factor, is a transcriptional activator that binds the CLCNKA promoter. KLF15, which is a nuclear protein, is expressed primarily in liver, heart, skeletal muscle and kidney tissues but is not detected in lymphoid tissues or bone marrow. It is an important regulator of GLUT4 in both adipose and muscle tissues.

抗原: Recombinant protein of human KLF15

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 µ g/ml.

特异性: Highly expressed in liver, followed by heart, skeletal muscle, and kidney. Not expressed in bone marrow or lymphoid tissues.