

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MDK Polyclonal Antibody

多克隆抗体

产品货号: PA06981

第三版

描述: This gene encodes a member of a small family of secreted growth factors that binds heparin and responds to retinoic acid. The encoded protein promotes cell growth, migration, and angiogenesis, in particular during tumorigenesis. This gene has been targeted as a therapeutic for a variety of different disorders. Alternatively spliced transcript variants encoding multiple isoforms have been observed.

抗原: Recombinant protein of human MDK

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in various tumor cell lines. In insulinoma tissue predominantly expressed in precancerous lesions.