

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human WDR36 Polyclonal Antibody

多克隆抗体

产品货号: PA07139

第三版

描述:This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asn (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Mutations in this gene have been associated with adult-onset primary open-angle glaucoma (POAG).

抗原:Synthetic peptide of human WDR36

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart, placenta, liver, skeletal muscle, kidney and pancreas. In ocular tissues, strong expression in iris, sclera, ciliary muscle, ciliary body, retina and optic nerve. 1 Publication