

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KDM5C Polyclonal Antibody

多克隆抗体

产品货号: PA07191

第三版

描述: This gene is a member of the SMCY homolog family and encodes a protein with one ARID domain, one JmjC domain, one JmjN domain and two PHD-type zinc fingers. The DNA-binding motifs suggest this protein is involved in the regulation of transcription and chromatin remodeling. Mutations in this gene have been associated with X-linked mental retardation. Alternative splicing results in multiple transcript variants.

抗原: Synthetic peptide of human KDM5C

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in all tissues examined. Highest levels found in brain and skeletal muscle.