

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HAUS7 Polyclonal Antibody

多克隆抗体

产品货号: PA07202

第三版

**描述:** This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants.

**抗原:** Synthetic peptide of human HAUS7

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Detected in spleen, thymus, testis, ovary, small intestine and colon, with highest levels of expression in testis and ovary.