

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TUSC1 Polyclonal Antibody

多克隆抗体

产品货号: PA07221

第三版

描述: This gene is located within the region of chromosome 9p that harbors tumor suppressor genes critical in carcinogenesis. It is an intronless gene which is downregulated in non-small-cell lung cancer and small-cell lung cancer cell lines, suggesting that it may play a role in lung tumorigenesis.

抗原: Synthetic peptide of human TUSC1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: