

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TSPY1/3 Polyclonal Antibody

多克隆抗体

产品货号: PA07232

第三版

描述:The protein encoded by this gene is found only in testicular tissue and may be involved in spermatogenesis. Approximately 35 copies of this gene are present in humans, but only a single, nonfunctional orthologous gene is found in mouse. Two transcript variants encoding different isoforms have been found for this gene.

抗原:Synthetic peptide of human TSPY1/TSPY3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: