
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TRIB1 Polyclonal Antibody

多克隆抗体

产品货号: PA07242

第三版

描述: Interacts with Human, Mouse APK kinases and regulates activation of MAP kinases. May not display kinase activity.

抗原: Synthetic peptide of human TRIB1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expressed in most human tissues with the highest levels in skeletal muscle, thyroid gland, pancreas, peripheral blood leukocytes, and bone marrow.