

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TRIM71 Polyclonal Antibody

多克隆抗体

产品货号: PA07247

第三版

描述: LIN-41, also called tripartite motif-containing 71 (TRIM71), which was first identified in *Caenorhabditis elegans*, is responsible for the timing of cell fate determination. By encoding microRNAs (miRNAs), the heterochronic genes *let-7* and *lin-4* downregulate the gene encoding LIN-41. The miRNAs bind to six complementary sites on the 3' untranslated region (UTR) of the LIN-41 gene. This downregulation positively regulates the timing of the expression of LIN-29, an adult specification transcription factor.

抗原: Synthetic peptide of human TRIM71

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Specifically expressed in testis.