

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC6A19 Polyclonal Antibody

多克隆抗体

产品货号: PA07270

第三版

描述: This gene encodes a system B(0) transmembrane protein that actively transports most neutral amino acids across the apical membrane of epithelial cells. Mutations in this gene result in Hartnup disorder.

抗原: Synthetic peptide of human SLC6A19

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Robust expression in kidney and small intestine, with minimal expression in pancreas. Also expressed in stomach, liver, duodenum, ileocecum, colon and prostate. Not detected in testis, whole brain, cerebellum, fetal liver, spleen, skeletal muscle, uterus, heart or lung.