

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DEFB112 Polyclonal Antibody

多克隆抗体

产品货号: PA07277

第三版

**描述:**Defensins form a family of antimicrobial and cytotoxic peptides made by neutrophils. Defensins are short, processed peptide molecules that are classified by structure into three groups: alpha-defensins, beta-defensins and theta-defensins. All beta-defensin genes are densely clustered in four to five syntenic chromosomal regions.

**抗原:**Synthetic peptide of human DEFB112

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**