



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NUF2 Polyclonal Antibody

多克隆抗体

产品货号: PA07306

第三版

描述: This gene encodes a protein that is highly similar to yeast Nuf2, a component of a conserved protein complex associated with the centromere. Yeast Nuf2 disappears from the centromere during meiotic prophase when centromeres lose their connection to the spindle pole body, and plays a regulatory role in chromosome segregation.

抗原: Synthetic peptide of human NUF2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: