

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EDG5 Polyclonal Antibody

多克隆抗体

产品货号: PA07311

第三版

描述: This gene encodes a member of the G protein-coupled receptors, as well as the EDG family of proteins. This protein participates in sphingosine 1-phosphate-induced cell proliferation, survival, and transcriptional activation

抗原: Synthetic peptide of human S1PR2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: