

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse USP39 Polyclonal Antibody

多克隆抗体

产品货号: PA07466

第三版

**描述:** USP39 (ubiquitin specific peptidase 39), also known as SAD1, CGI-21, HSPC332 or SNRNP65, is a 565 amino acid nuclear protein belonging to the peptidase C19 family. Containing one UBP-type zinc finger, USP39 play a role in mRNA splicing and is part of a tri-snRNP complex. USP39 is thought to be required for the stabilization of the spindle checkpoint and may help support successful cytokinesis.

**抗原:** Fusion protein of human USP39

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:**