

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human UGT1A10 Polyclonal Antibody

多克隆抗体

产品货号: PA07492

第三版

描述: This gene encodes a UDP-glucuronosyltransferase, an enzyme of the glucuronidation pathway that transforms small lipophilic molecules, such as steroids, bilirubin, hormones, and drugs, into water-soluble, excretable metabolites. This gene is part of a complex locus that encodes several UDP-glucuronosyltransferases. The locus includes thirteen unique alternate first exons followed by four common exons. Four of the alternate first exons are considered pseudogenes.

抗原: Fusion protein of human UGT1A10

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Liver and colon. Isoform 1 and isoform 2 are expressed in colon, esophagus and small intestine; isoform 2 but not isoform 1 is expressed in liver or kidney (PubMed:18004212).