

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NME8 Polyclonal Antibody

多克隆抗体

产品货号: PA07516

第三版

**描述:** This gene encodes a protein with an N-terminal thioredoxin domain and three C-terminal nucleoside diphosphate kinase (NDK) domains, but the NDK domains are thought to be catalytically inactive. The sea urchin ortholog of this gene encodes a component of sperm outer dynein arms, and the protein is implicated in ciliary function. Mutations in this gene are implicated in primary ciliary dyskinesia type 6.

**抗原:** Fusion protein of human NME8

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Testis-specific. Expressed only in primary spermatocytes and round spermatids. 1 Publication