

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NME8 Polyclonal Antibody

多克隆抗体

产品货号: PA07517

第三版

描述: This gene encodes a protein with an N-terminal thioredoxin domain and three C-terminal nucleoside diphosphate kinase (NDK) domains, but the NDK domains are thought to be catalytically inactive. The sea urchin ortholog of this gene encodes a component of sperm outer dynein arms, and the protein is implicated in ciliary function. Mutations in this gene are implicated in primary ciliary dyskinesia type 6.

抗原: Fusion protein of human NME8

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Testis-specific. Expressed only in primary spermatocytes and round spermatids. 1 Publication