

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TULP2 Polyclonal Antibody

多克隆抗体

产品货号: PA07522

第三版

**描述:**TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues.

**抗原:**Fusion protein of human TULP2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Strongly expressed in testis. Also expressed in retina. Expressed in cancer cell lines.