



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRAK1 Polyclonal Antibody

多克隆抗体

产品货号: PA07558

第三版

描述:Involved in the regulation of endosome-to-lysosome trafficking, including endocytic trafficking of EGF-EGFR complexes and GABA-A receptors.

抗原:Fusion protein of human TRAK1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:High expression in spinal cord and moderate expression in all other tissues and specific brain regions examined. Expressed in all cell lines examined.1 Publication