

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRAK1 Polyclonal Antibody

多克隆抗体

产品货号: PA07559

第三版

**描述:**Involved in the regulation of endosome-to-lysosome trafficking, including endocytic trafficking of EGF-EGFR complexes and GABA-A receptors.

**抗原:**Fusion protein of human TRAK1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**High expression in spinal cord and moderate expression in all other tissues and specific brain regions examined. Expressed in all cell lines examined.1 Publication