

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TMPRSS11D Polyclonal Antibody

多克隆抗体

产品货号: PA07569

第三版

描述: This gene encodes a trypsin-like serine protease released from the submucosal serous glands onto mucous membrane. It is a type II integral membrane protein and has 29-38% identity in the sequence of the catalytic region with human hepsin, enteropeptidase, acrosin, and mast cell tryptase. The noncatalytic region has little similarity to other known proteins.

抗原: Fusion protein of human TMPRSS11D

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Located in the cells of the submucosal serous glands of the bronchi and trachea.