
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TMPRSS11E Polyclonal Antibody

多克隆抗体

产品货号: PA07572

第三版

描述: Serine protease which possesses both gelatinolytic and caseinolytic activities. Shows a preference for Arg in the P1 position (By similarity).

抗原: Fusion protein of human TMPRSS11E

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expression can only be detected in tissues derived from the head and neck, and in skin, prostate and testis.