

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TMPRSS11E Polyclonal Antibody

多克隆抗体

产品货号: PA07573

第三版

**描述:** Serine protease which possesses both gelatinolytic and caseinolytic activities. Shows a preference for Arg in the P1 position (By similarity).

**抗原:** Fusion protein of human TMPRSS11E

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expression can only be detected in tissues derived from the head and neck, and in skin, prostate and testis.