

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TMPRSS13 Polyclonal Antibody

多克隆抗体

产品货号: PA07575

第三版

描述: This gene encodes a member of the type II transmembrane serine protease family. Transmembrane serine proteases are regulated by protease inhibitors and known to function in development, homeostasis, infection, and tumorigenesis. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原: Fusion protein of human TMPRSS13

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Isoform 1 and isoform 3 are predominantly expressed in lung, placenta, pancreas, and prostate. Isoform 3 is weakly expressed in testis and peripheral blood lymphocytes.