

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse Tmprss5 Polyclonal Antibody

多克隆抗体

产品货号: PA07576

第三版

**描述:** This gene encodes a protein that belongs to the serine protease family. Serine proteases are known to be involved in many physiological and pathological processes. Alternative splicing results in multiple transcript variants.

**抗原:** Fusion protein of human Tmprss5

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Brain-specific. Predominantly expressed in neurons, in their axons, and at the synapses of motoneurons in the spinal cord.