

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse RASD2 Polyclonal Antibody

多克隆抗体

产品货号: PA07626

---

第三版

**描述:** This gene encodes a Ras-related protein that enriched in striatum. The product of this gene binds to GTP and possesses intrinsic GTPase activity. The gene belongs to the Ras superfamily of small GTPases. The exact function of this gene is unknown, but most striatum-specific mRNAs characterized to date encode components of signal transduction cascades.

**抗原:** Full length fusion protein

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Pancreatic endocrine cells (islets of Langerhans).