

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat STX12 Polyclonal Antibody

多克隆抗体

产品货号: PA07639

第三版

**描述:** Syntaxin 12, also known as STX12, STX13 or STX14, is a 276 amino acid single-pass membrane protein that contains one t-SNARE coiled-coil homology domain and belongs to the Syntaxin family. Syntaxin 12 regulates protein transport between late endosomes and the trans-Golgi network and interacts with ABC1 (ATP-binding cassette transporter A1), a protein that facilitates cellular release of choline-phospholipids and cholesterol to apoA-I (apolipoprotein A-I).

**抗原:** Full length fusion protein

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:**