

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ECE2 Polyclonal Antibody

多克隆抗体

产品货号: PA07699

第三版

描述: This gene encodes a member of the M13 family, which includes type 2 integral membrane metalloproteinases. The encoded enzyme is a membrane-bound zinc-dependent metalloprotease. The enzyme catalyzes the cleavage of big endothelin to produce the vasoconstrictor endothelin-1, and plays a role in the processing of several neuroendocrine peptides. It may also have methyltransferase activity. Alternative splicing results in multiple transcript variants.

抗原: Fusion protein of human ECE2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in brain. Strongly down-regulated in inferior parietal lobe from Alzheimer disease patients (at protein level).