

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GLP2R Polyclonal Antibody

多克隆抗体

产品货号: PA07715

第三版

描述: This gene encodes a G protein-coupled receptor that is closely related to the glucagon receptor and binds to glucagon-like peptide-2 (GLP2). Signalling through GLP2 stimulates intestinal growth and increases villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis.

抗原: Fusion protein of human GLP2R

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: