

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse GLP2R Polyclonal Antibody

多克隆抗体

产品货号: PA07716

---

第三版

**描述:** This gene encodes a G protein-coupled receptor that is closely related to the glucagon receptor and binds to glucagon-like peptide-2 (GLP2). Signalling through GLP2 stimulates intestinal growth and increases villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis.

**抗原:** Fusion protein of human GLP2R

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**