

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRPM6 Polyclonal Antibody

多克隆抗体

产品货号: PA07733

第三版

描述: This gene is predominantly expressed in the kidney and colon, and encodes a protein containing an ion channel domain and a protein kinase domain. It is crucial for magnesium homeostasis, and plays an essential role in epithelial magnesium transport and in the active magnesium absorption in the gut and kidney. Mutations in this gene are associated with hypomagnesemia with secondary hypocalcemia. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

抗原: Synthetic peptide of human TRPM6

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Highly expressed in kidney and colon. Isoform TRPM6a and isoform TRPM6b, are coexpressed with TRPM7 in kidney, and testis, and are also found in several cell lines of lung origin. Isoform TRPM6c is detected only in testis and in NCI-H510A small cell lung carcinoma cells.