

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PKA gamma Polyclonal Antibody

多克隆抗体

产品货号: PA07756

第三版

描述: cAMP-dependent protein kinase catalytic subunit gamma is an enzyme that in humans is encoded by the PRKACG gene. Cyclic AMP-dependent protein kinase (PKA) consists of two catalytic subunits and a regulatory subunit dimer. This gene encodes the gamma form of its catalytic subunit. The gene is intronless and is thought to be a retrotransposon derived from the gene for the alpha form of the PKA catalytic subunit. PRKACG has been shown to interact with Ryanodine receptor 2.

抗原: Recombinant protein of human PRKACG

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Testis specific. But important tissues such as brain and ovary have not been analyzed for the content of transcript.