

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LGALS7 Polyclonal Antibody

多克隆抗体

产品货号: PA07784

第三版

描述:The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. Differential and in situ hybridization studies indicate that this lectin is specifically expressed in keratinocytes and found mainly in stratified squamous epithelium. A duplicate copy of this gene (GeneID:653499) is found adjacent to, but on the opposite strand on chromosome 19.

抗原:Synthetic peptide of human LGALS7

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Mainly expressed in stratified squamous epithelium.