

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CMTM7 Polyclonal Antibody

多克隆抗体

产品货号: PA07791

第三版

**描述:** This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. The protein encoded by this gene is highly expressed in leukocytes, but its exact function is unknown. Alternatively spliced transcript variants encoding different isoforms have been identified.

**抗原:** Synthetic peptide of human CMTM7

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Highly expressed in leukocytes.1 Publication