

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CMTM6 Polyclonal Antibody

多克隆抗体

产品货号: PA07793

第三版

**描述:**This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. This gene is widely expressed in many tissues, but the exact function of the encoded protein is unknown.

**抗原:**Synthetic peptide of human CMTM6

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in the leukocytes, placenta and testis.