

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SYT9 Polyclonal Antibody

多克隆抗体

产品货号: PA07829

第三版

描述: Synaptotagmin-9 is a protein that in humans is encoded by the SYT9 gene. SYT9 has been shown to interact with SYNCRIP, TUBB and TRPV1. May be involved in Ca²⁺-dependent exocytosis of secretory vesicles through Ca²⁺ and phospholipid binding to the C2 domain or may serve as Ca²⁺ sensors in the process of vesicular trafficking and exocytosis.

抗原: Recombinant protein of human SYT9

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: